

PASS Labs X250.5

Masuratori la data INITIAL : *cu capacul deschis în echilibru termic*

LEFT ch.						RIGHT ch.					
TRANZ. model	REZ. 0,47 Ohm			REZ. 221 Ohm		TRANZ. model	REZ. 0,47 Ohm			REZ. 221 Ohm	
	Nr.	Ohm	mV	Ohm	mV		Nr.	Ohm	mV	Ohm	mV
IRF9240	R38	0.4	7,7	221		IRF240	R58	0.4	49,5	221	
IRF9240	R30	0.4	6,1	222		IRF240	R59	0.4	48	222	
IRF9240	R22	0.4	6,6	221		IRF240	R61	0.4	46,5	222	
IRF240	R37	0.4	47,5	220		IRF9240	R63	0.4	9,5	221	
IRF240	R29	0.4	47,5	219		IRF9240	R65	0.4	8,1	222	
IRF240	R21	0.4	44,5	219		IRF9240	R67	0.4	8	222	
IRF9240	R36	0.4	20,5	219		IRF240	R85	0.4	43,5	221	
IRF9240	R28	0.4	20	219		IRF240	R87	0.4	42,3	221	
IRF9240	R20	0.4	20,5	219		IRF240	R71	0.4	42	222	
IRF240	R35	0.4	59,5	222		IRF9240	R72	0.4	5,7	220	
IRF240	R27	0.4	60	221		IRF9240	R74	0.4	0,4	219	
IRF240	R19	0.4	60 (62)	221		IRF9240	R76	0.4	5,5	219	

Temperatura de echilibru 39,1 °C Temperatura de echilibru 38 °C

